

In-Chassis Calibration of PMMs

16 November 2020

M. Daniels

This tutorial addresses how to automate the calibration and adjustment of the Fluke Calibration Pressure Measurement Modules (“PMM”) while they remain in the chassis of a 6270A, 8270A, or 8370A. The automation scheme is similar to that of the PMM calibration when used with the Cal Sled accessory – (link: <https://support.flukecal.com/hc/en-us/articles/115005079666-Macros-PM200-and-PM600-Fully-Automated-Calibration>).

Like the Cal Sled concept, the In-Chassis calibration provides a fully automated process beginning with identifying the PMM to be calibrated, collecting As-Found data, making a correction based on the A/F data, and then completing the process with an As-Left verification data set. There is an optional guard banding assessment which can be used to determine final acceptance of the results. The In-Chassis calibration is applicable to all PMMs regardless of type (PM200 / PM500 / PM600) or pressure range.

The automation is attained using COMPASS for Pressure macros. They were built and tested using COMPASS for Pressure v5.0.50. An “Enhanced” license for the software is required.



Table of Contents: (Click the embedded link to jump to that section of the tutorial)

- I. [Macro Overview](#)
- II. [COMPASS Device Definition Setups](#)
- III. [COMPASS Options](#)
- IV. [Initializing the Calibration](#)
- V. [Performing the Calibration](#)
- VI. [Restoring the Previous Calibration](#)
- VII. [Test Complete](#)
- VIII. [Macro Code](#)

I. Macro Overview

There are four Test Macros required for the In-Chassis Calibrations. The macros can be manually added or imported by using the provided links at the end of this tutorial.

- ChassisCal_Master
- ChassisCal_PreTest
- ChassisCal_Absolute
- ChassisCal_GaugeDiff

“ChassisCal_Master”: this macro is assigned as a Test Event macro as part of the COMPASS Test Definition. It is used to call the main automation macro based on the measurement mode of the PMM. It also contains subroutines used to set the calibration date and to restore the calibration should the As Left results not be desired.

“ChassisCal_PreTest”: this macro is assigned as a Pre-Test macro as part of the COMPASS Test Definition. It queries the chassis for which PMMs are installed. Then it prompts the operator to select the one to be calibrated. The slot address and slot number ID are assigned and used throughout the rest of the calibration to build device-specific commands. The PMM is fixed as the active module in the chassis and identified in the data file as the one being calibrated. This enables the data directory (ie: “C:\dhi\COMPASS for Pressure\Data...”) to be named by the serial number of the PPM being calibrated and not by the chassis serial number. Lastly, the fluid head is set to zero.


“ChassisCal_Absolute”: this macro is not assigned as part of a Test Definition, but rather it is called by the ChassisCal_Master and runs in the background. It is the main macro used to effect the adjustment of the PMMs operating in absolute mode. For example, calibrating a PM600-A3.5M.

“ChassisCal_GaugeDiff”: this macro is not assigned as part of a Test Definition, but rather it is called by the ChassisCal_Master and runs in the background. It is the main macro used to effect the adjustment of the PMMs operating in gauge or differential mode. For example, calibrating a PM200-G2M, or a PM500-BG1M.


II. COMPASS Device Definition Setups

COMPASS Device Definitions, DUT:

The DUT Definitions will follow the normal configuration methods of an Auto-Detect device. It is anticipated that most users will already have the 6270A or 8270A / 8370A configured as a Support Device. If so, then the checkbox at the bottom of the Header tab can be enabled to make the device definition useable as a DUT.



The Tolerance tab will be used by COMPASS to determine the In or Out of Tolerance status of the PMM being calibrated. The data used on the Tolerance tab needs to reflect the tolerance of the PMM being calibrated. In the situation where different types of PMMs are used in the same chassis then a decision needs to be made by each organization on how the Tolerancing will be handled. One option is to edit this information as part of the *Customize Output* during the test initialization sequence. Another option is to create individual Device Definitions for each PMM, and then have the Record Label indicate a unique identification. This tutorial will demonstrate the *Customize Output* method during test initialization. The screen capture below is for a PM600 module.




COMPASS Device Definitions, Support Device:


The Support Device definition can be any device deemed appropriate for the calibration of a PMM. Normal methods of configuring the Support Device definition apply.

COMPASS Device Definitions, Test:


The Test Definition is the key to making the automated process work. This is where the macros are assigned and the two pressure cycles enabled. The Pre Test Macro “ChassisCal_PreTest” must be assigned on the Pre-Test tab.



The Pressure tab requires 2 Pressure Cycles be assigned. The first cycle is for the As Found data, and the As Left data is the second cycle. The calibration coefficients are modified between the two cycles. A typical calibration sequence may look like the following...




The Data tab is where the ChassisCal_Master macro must be assigned as a Test Event Macro:




No Auxiliary devices are required for automated calibration.

III. COMPASS Options

There are three selections required from the [Tools],<Options> menu. In the “Data in File” tab the fields named “String Data 1”, “String Data 2” and “String Data 3” are required to be in the *Fields to log* category. These fields are used to store the PMM slot ID, Address and type.




In the “Initialize Test” tab, the option to *Load calibration info for auto-detect devices* must be disabled by unchecking the box.



IV. Initializing the Calibration:

In these screen shots, a PM200-G2M is being calibrated and is located in a 8270A chassis. The reference device is a PM600 and is located in a 6270A chassis. This arrangement is made for the convenience of demonstrating the process and is not an endorsement of metrological practices in terms of Test Accuracy Ratios.




When the instrument capabilities have been uploaded to COMPASS, then the operator is required to, 1) define the range of the PMM being calibrated, and 2) customize the Tolerance information to that of the PMM being calibrated.

- Define the range of the PPM to be calibrated: A key part of the In-Chassis calibration is the ability to lock in on the PMM to be calibrated. To this end, the PreTest macro will prompt the operator to (again) select the PMM to be calibrated. It's the macro-based selection which locks in the range and fixes the PMM to be used for the calibration. However, the Min and Max fields are important if the Test Definition is defined in units "% DUT Span". The option to select "Auto Range" is a valid choice.




For gauge mode and differential mode calibrations it is important to check the box for "Enable Auto Zero". This is the mechanism which will allow the chassis to zero ("tare") the PMM when vented. For absolute mode PMMs the implications are less significant and depends on the presence and utilization of a barometric range module. The checkbox acts as an On/Off switch for the chassis to auto zero or not to auto zero. Regardless of the state of the checkbox the measured pressure will always use the PMM's zOffset coefficient.


- Customize the Tolerance information: Although a standard feature, this method is not commonly used in COMPASS. It provides a great deal of flexibility for defining device outputs. The tolerance information is changed by clicking the blue-font text "Customize Output". This selection is important for COMPASS to determine In or Out of Tolerance, and for the correct display of the plot tolerance bars. For the PM200-G2M the tolerance is $\pm 0.02\%$ FS.




A Test Definition must be selected...




Select the Reference and Control devices...



Set the Fluid Head height if necessary...




Final confirmation of the test to be performed... Click Finish to begin the calibration.




Performing the Calibration:

When the initialization steps have been completed then the automated test begins. The PreTest macro, "ChassisCal_PreTest" scans the modules installed in the chassis and prompts the operator to select the PMM to be calibrated. This user selection is used by COMPASS to make the desired PMM the "fixed" selection.



After the PMM selection is made then COMPASS progresses through the steps of the Test Definition. At the end of the first pressure cycle the data is automatically evaluated using a linear regression and new coefficients are determined and written to the PMM. The As Left verification cycle is then automatically started.




V. Restore As-Found coefficients:

The last part of the automated testing is an evaluation of the As-Left data to determine if the corrected results are inside a user-defined acceptance tolerance. The default criterion is 70% of the PMM standard specification. This is a form of guardbanding. If a different acceptance criterion is desired then the “TOLFACT” variable near the end of the ChassisCal_GaugeDiff or ChassisCal_Absolute macros can be edited. In the image below a TOLFACT of .70 is used.


```
10193 ****
10194
10195 'This block of code is used to check the fit of the As Left corrections. The key variable is named "TOLFACT".
10196 'If TOLFACT = 1 then the As Left data must fit inside the 1 year tolerance specification. If set to <1 then
10197 'this acts as a form of guardbanding.
10198
10199 TOLFACT = .70    '70% guardband. User definable.
10200
10201 Dim MessageBox1
10202 Dim MessageBox2
10203 For i = 1 To cCOMPASS.DataCollection(1).NumberofPressurePoints
```

If any part of the data set does not fall inside of the 70% guardband then a messagebox is generated alerting the user to this condition. The option is given to revert (restore) to the original As Found coefficients. The default option is to keep the newly determined coefficients. If the choice is made to restore then a second prompt is given to make sure this is the desired action. Restoring will undo the results of the test.



VI. Test Complete


When the automated testing is completed then a post test dialog box will appear. Following standard COMPASS conventions, the data file is always saved after each test point is logged.



VII. Macro Code


There are two methods for placing the macros into [your] active database – import, or copy-and-paste.

- Import method: the macros and example device definitions are provided as a link at the bottom of this tutorial. The link downloads a COMPASS database “compsetp_5.1.mdb”. The content can be transferred in your active database using the COMPASS Database Maintenance Tool.



- b) Copy-and-Paste method: VB script code is provided below which can be copied and then pasted into the COMPASS Macro Editor. New Macro folders will be required before the code can be copied.

The macros are Test macros.



The screenshot shows the COMPASS Macro Editor interface. The title bar says "COMPASS Macro Editor". The menu bar has "Edit" and "Settings". The toolbar includes icons for file operations and help. The left pane shows a tree view of macro files under "All Code". The "Test" folder is expanded, containing sub-macros like StepBack, WriteTest, 34401NULL, 7252LoopTest, AbsMode, AutoZeroRPM4AD, BGNegativeExercise, CalActivate, CalActivate_PreTest, CalibrateADC2500, CalibrateADC2500PrePoint, ChassisCal_Absolute, ChassisCal_GaugeDiff, ChassisCal_Master (which is selected), ChassisCal_PreTest, CloseVent, COMPASS formatted MessageBox, CustomerInformation, CustomerInformation2, DeactivateValveDriver, DHallMacroExample, EDWTCalEvent, Elmendorf_Macro, EventSound, Fluke700_700P_CalibrationAdjust_DV, and Fluke700_700P_CalibrationAdjust_DV. The right pane shows the code for the selected macro, ChassisCal_Master:

```

Title: ChassisCal_Master

9799 |
9800 'Used with the 6270A / 8x70A in-chassis calibration of the PMMs.
9801 'The main TestEvent Macro to call. It's used to call one of two
9802 'separate Test Macros based on the measurement mode of the PMM.
9803 ****
9804
9805 Function ChassisCalibrationAdjust(iT, iL, iC, iP, cTest, cConfig)
9806
9807 'Determine which macro to run. There are five Meas Mode codes:
9808 ' 0 = Absolute
9809 ' 1 = Gauge
9810 ' 2 = Abs by Atms
9811 ' 3 = Differential (Bi-Directional)
9812 ' 4 = ???
9813 ' 5 = Negative Gauge
9814
9815 Select Case cTest.TestPrsMeasMode
9816     Case 0:
9817         Call ChassisCal_Absolute(iT, iL, iC, iP, cTest, cConfig)
9818     Case 1:
9819         Call ChassisCal_GaugeDiff(iT, iL, iC, iP, cTest, cConfig)
9820     Case 3:
9821         Call ChassisCal_GaugeDiff(iT, iL, iC, iP, cTest, cConfig)

```

ChassisCal_Master:

```
'Used with the 6270A / 8x70A in-chassis calibration of the PMMs.  
'The main TestEvent Macro to call. It's used to call one of two  
'separate Test Macros based on the measurement mode of the PMM.  
*****  
Function ChassisCalibrationAdjust(iT, iL, iC, iP, cTest, cConfig)  
  
'Determine which macro to run. There are five Meas Mode codes:  
'    0 = Absolute  
'    1 = Gauge  
'    2 = Abs by Atms  
'    3 = Differential (Bi-Directional)  
'    4 = ???  
'    5 = Negative Gauge  
  
    Select Case cTest.TestPrsMeasMode  
        Case 0:  
            Call ChassisCal_Absolute(iT, iL, iC, iP, cTest, cConfig)  
        Case 1:  
            Call ChassisCal_GaugeDiff(iT, iL, iC, iP, cTest, cConfig)  
        Case 3:  
            Call ChassisCal_GaugeDiff(iT, iL, iC, iP, cTest, cConfig)  
  
    End Select  
'NOTE: Case 3 is the "Differential" selection in the Test Definition. The intent  
'is to start both gauge and differential mode calibrations at a vented condition  
'and replicate the functionality of the user pressing the "zero" button. Once the  
'PMM is zeroed then the test moves to set the first target pressure.  
  
End Function  
  
Function ChassisCal_Date()  
    YY = right(Year(Date),2)  
    dd = Right(String(2, "0") & Day(Date), 2)  
    mm = Right(String(2, "0") & Month(Date), 2)  
    ChassisCal_Date = "20" & YY & "," & mm & "," & dd  
End Function  
  
'This is called by ChassisCal_Absolute and ChassisCal_GaugeDiff to restore  
'the original As Found coefficients.  
  
Function RestoreCal()  
cdebug.LogStatus "The RestoreCal function has been called."  
  
'Query the original A/F coefficients:  
origC0 = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient1  
origC1 = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient2  
origzOffset = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient3  
old_date = cCOMPASS.DataCollection(1).DUT.LastCalDate  
'Convert to six fixed decimal places:  
OriginalC0 = FormatNumber(origC0,6)  
OriginalC1 = FormatNumber(origC1,6)  
OriginalzOff = FormatNumber(origzOffset,6)  
    cdebug.LogStatus "the original C0 is: " & OriginalC0  
    cdebug.LogStatus "the original C1 is: " & OriginalC1  
    cdebug.LogStatus "the original zOffset is: " & OriginalzOff  
    cDebug.LogStatus "Calibration & Adjustment date: " & old_date  
  
'Query the new A/L coefficients:  
newC0 = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient4  
newC1 = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient5  
newzOffset = cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient6  
'Convert to six fixed decimal places:  
NewC0six = FormatNumber(newC0,6)  
NewC1six = FormatNumber(newC1,6)  
NewzOffsix = FormatNumber(newzOffset,6)  
    cdebug.LogStatus "the A/L 'NewC0six' coef is: " & NewC0six  
    cdebug.LogStatus "the A/L 'NewC1six' coef is: " & NewC1six  
    cdebug.LogStatus "the A/L 'NewzOffsix' coef is: " & NewzOffsix  
  
'Build query commands to help confirm the original A/F coefs have been properly written:  
PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2  
cmd1 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1?"  
cmd2 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2?"  
cmd3 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"  
cmd4 = "CAL:PRES" & PMMSlotAddress & ":DATE?"  
  
'Build write commands to restore the original A/F coefs:  
cmd5 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1" & OriginalC0  
cmd6 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2" & OriginalC1  
cmd7 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA" & OriginalzOff  
cmd8 = "CAL:PRES" & PMMSlotAddress & ":DATE" & old_date
```

```

cdebug.LogStatus "cmd5 is: " & cmd5
cdebug.LogStatus "cmd6 is: " & cmd6
cdebug.LogStatus "cmd7 is: " & cmd7
cdebug.LogStatus "cmd8 is: " & cmd8

Do
    'Enable the CAL Mode:
    cConfig.DUTPrs(1).IoSendCommand "CAL:MODE 1",True
    TimeDelay 1
    mode = cConfig.DUTPrs(1).IoSendCommand ("CAL:MODE?",False)
    cdebug.LogStatus "Cal Mode is: " & mode

    'Write to restore the original coeffs:
    Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd5),True)
    Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd6),True)
    Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd7),True)
    Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd8),True)

    'Read back the written coeffs:
    ReadBackC0raw = cConfig.DUTPrs(1).IoSendCommand (CStr(cmd1),False)
    ReadBackClraw = cConfig.DUTPrs(1).IoSendCommand (CStr(cmd2),False)
    ReadBackZ = cConfig.DUTPrs(1).IoSendCommand (CStr(cmd3),False)
        cdebug.LogStatus "ReadBackC0raw value is: " & ReadBackC0raw
        cdebug.LogStatus "ReadackClraw value is: " & ReadBackClraw
        cdebug.LogStatus "ReadBackZ value is: " & ReadBackZ
    'Remove the leading coef identifier to attain only the value...
    ReadBackC0 = qextract(ReadBackC0raw,1,2,"")
    ReadBackCl = qextract(ReadBackClraw,1,2,"")
        cDebug.LogStatus "ReadBackC0 value is: " & ReadBackC0
        cDebug.LogStatus "ReadBackCl value is: " & ReadBackCl

    'Check if the restored coef is different from the new A/L (it should be):
    'But first, limit each value to six decimal places...
    ReadbackC0six = FormatNumber(ReadBackC0,6)
    ReadbackClsix = FormatNumber(ReadBackCl,6)
    ReadBackZsix = FormatNumber(ReadBackZ,6)
        cdebug.LogStatus "ReadbackC0six value is: " & ReadBackC0six
        cdebug.LogStatus "ReadbackClsix value is: " & ReadBackClsix
        cdebug.LogStatus "ReadBackZsix value is: " & ReadBackZsix

    C0Diff = NewC0six - ReadBackC0six
    C1Diff = NewClsix - ReadBackClsix
    ZoffDiff = NewZoffsix - ReadBackZsix
        cDebug.LogStatus "C0 restore Difference: " & C0Diff
        cDebug.LogStatus "C1 restore Difference: " & C1diff
        cDebug.LogStatus "Zoff restore Difference: " & ZoffDiff

    If cCOMPASS.SystemAbort Then Exit Function
    If C0Diff <> 0 And C1Diff <> 0 Then Exit Do
        cdebug.LogStatus "End of Loop."

'NOTE: PM500 modules have the C0 reported in "kPa" and displayed on the chassis
'front panel in units of "psi" in the latest versions of chassis firmware.
'vel.10 = 6720A
'vel.02 = 8x70A

cCOMPASS.StatusDisplay "Restoring the original As Found coefficients..."
Loop

End Function

```

ChassisCal_PreTest:

```

'*****
'This Pre-Test macro queries the chassis for which PMMs are installed.
'Then it prompts the operator to select the one to be calibrated. The PMM is
'fixed as the active one, identified in the .dat file as the one being calibrated,
'and the fluid head set to 0.
'*****
Function GatherModuleInfo(iT, iL, iC, iP, cTest, cConfig)

'scan each slot and return the ID and serial number:
PMM1 = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD1:NAME?",False,0)
PMM1sn = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD1:SER?",False,0)
    cDebug.LogStatus "PMM1 Name: " & PMM1
    cDebug.LogStatus "PMM1 S/N: " & PMM1sn
PMM2 = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD2:NAME?",False,0)
PMM2sn = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD2:SER?",False,0)
    cDebug.LogStatus "PMM2 Name: " & PMM2
    cDebug.LogStatus "PMM2 S/N: " & PMM2sn
PMM3 = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD3:NAME?",False,0)
PMM3sn = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD3:SER?",False,0)
    cDebug.LogStatus "PMM3 Name: " & PMM3
    cDebug.LogStatus "PMM3 S/N: " & PMM3sn
PMM4 = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD4:NAME?",False,0)
PMM4sn = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD4:SER?",False,0)
    cDebug.LogStatus "PMM4 Name: " & PMM4
    cDebug.LogStatus "PMM4 S/N: " & PMM4sn
PMM5 = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD5:NAME?",False,0)
PMM5sn = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD5:SER?",False,0)
    cDebug.LogStatus "PMM5 Name: " & PMM5
    cDebug.LogStatus "PMM5 S/N: " & PMM5sn

'Present Inputbox message listing available PMMs, asking the user to make a selection:
With cConfig.DUTPrs(1).RangeMain.GetParent
    msg = "Enter the number of the PMM to be calibrated: " & vbCrLf
    msg = msg & vbCrLf
    msg = msg & "1. Model: " & PMM1 & vbCrLf
    msg = msg & "    s/n: " & PMM1sn & vbCrLf
    msg = msg & "2. Model: " & PMM2 & vbCrLf
    msg = msg & "    s/n: " & PMM2sn & vbCrLf
    msg = msg & "3. Model: " & PMM3 & vbCrLf
    msg = msg & "    s/n: " & PMM3sn & vbCrLf
    msg = msg & "4. Model: " & PMM4 & vbCrLf
    msg = msg & "    s/n: " & PMM4sn & vbCrLf
    msg = msg & "5. Model: " & PMM5 & vbCrLf
    msg = msg & "    s/n: " & PMM5sn

    reply = inputbox(msg,"PMMs Detected in Chassis", "")
    cDebug.LogStatus "The selection entered: " & reply

    'make sure a PMM was selected by the user:
    If reply = "" Then
        cdebug.LogStatus "A user input for the PMM selection was not detected. Aborted the test."
        msg = msgbox("A user input for the PMM selection was not detected. Aborting the test.", ,0+16+4096, "Test Aborted")
        cCOMPASS.SystemAbort = True
    End If

End With

'Based on the user input, select the PMM and fix it as active PMM:
If reply = 1 Then
    PMMSlotAddress = 5
    FixPMM = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD 1",True,0)
    cCOMPASS.TimeDelay 2
ElseIf reply = 2 Then
    PMMSlotAddress = 15
    FixPMM = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD 2",True,0)
    cCOMPASS.TimeDelay 2
ElseIf reply = 3 Then
    PMMSlotAddress = 25
    FixPMM = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD 3",True,0)
    cCOMPASS.TimeDelay 2
ElseIf reply = 4 Then
    PMMSlotAddress = 35
    FixPMM = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD 4",True,0)
    cCOMPASS.TimeDelay 2
ElseIf reply = 5 Then
    PMMSlotAddress = 45
    FixPMM = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD 5",True,0)
    cCOMPASS.TimeDelay 2
End If

```

```

PMMSlotNum = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD?",False,0)
cDebug.LogStatus "The fixed PMM Slot is: " & PMMSlotNum
cDebug.LogStatus "The address of the PMM to be calibrated: " & PMMSlotAddress

'SlotNumber = 1-5. SlotAddress = 5-45.
'write the slot and address of the PMM to the .dat file:
cConfig.DUTPrs(1).RangeMain.GetParent.Extra1 = PMMSlotNum
cConfig.DUTPrs(1).RangeMain.GetParent.Extra2 = PMMSlotAddress

'Set fluid head to 0. (COMPASS should have already done this as part of DUT initialization):
cConfig.DUTPrs(1).IoSendCommand "SENS:REF:HEIG 0",True,0

'Save the the .dat file which PMM was selected:
'Re-query the info:
Model = cConfig.DUTPrs(1).IoSendCommand ("SENS:MOD" & reply & ":NAME?",False)
SN = cconfig.DUTPrs(1).IoSendCommand ("SENS:MOD" & reply & ":SER?",False)
cdebug.LogStatus "the Model selected: " & Model
cdebug.LogStatus "the sn selected: " & SN

'write the info to the .dat file: (Note - this will be used by Windows to create / select
'a new sub directory based on the PMM's serial number)
cConfig.DUTPrs(1).RangeMain.GetParent.Model = Model
cConfig.DUTPrs(1).RangeMain.GetParent.SN = SN

'Define a Model Type variable for use in the unit conversion needs:
ModelType = Left(Model,5)
cdebug.LogStatus "the Model type selected is:" & ModelType

If ModelType = "PM200" Then
    ModelTypeID = 2
    cConfig.DUTPrs(1).RangeMain.GetParent.Extra3 = ModelTypeID
ElseIf ModelType = "PM500" Then
    ModelTypeID = 5
    cConfig.DUTPrs(1).RangeMain.GetParent.Extra3 = ModelTypeID
ElseIf ModelType = "PM600" Then
    ModelTypeID = 6
    cConfig.DUTPrs(1).RangeMain.GetParent.Extra3 = ModelTypeID
End If
cdebug.LogStatus "The ModelTypeID is: " & ModelTypeID
End Function

```

ChassisCal_GaugeDiff:

```
'*****
'Used with ChassisCal_Master to calibrate and adjust Gauge and Differential mode PMMs.
'Applicable for 6270A / 8270A / 8370A. Called by the master test macro "ChassisCal_Master".
'Gauge Mode testing must be performed with the user-option "Load calibration info for auto-detect
'devices" disabled.
'*****  
  
Function ChassisCal_GaugeDiff(iT, iL, iC, iP, cTest, cConfig)
Dim DUT_Raw_Data()
Dim Reference_Data()
Dim C1_Corrected_Data()
Dim res(4,7)  
  
C0 = 0
C1 = 0
zOffset = 0  
  
'<currentteststep 1010 - record AF coefs.>
If cCOMPASS.CurrentTestStep = 1010 Then ' Test Definition data files were created for each DUT.
    cDebug.LogStatus "Starting <step 1010> "
    cDebug.LogStatus "Capture As Found Coeffs..."  
  
'Ensure the Chassis is in Gauge mode:
cConfig.DUTPrs(1).IoSendCommand "SENS:MODE GAUGE",True  
  
PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2
cDebug.LogStatus "The PMM slot address is: " & PMMSlotAddress  
  
cmd1 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1?"
cmd2 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2?"
cmd3 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
cmd4 = "CAL:PRES" & PMMSlotAddress & ":DATE?"  
  
'Identify the unit of measure being used for the test and the PMM Type:
DUTUnit = cConfig.DUTPrs(1).RangeMain.UnitFinal
DUTunitText = cConfig.DUTPrs(1).RangeMain.UnitFinalText
ModelTypeID = cConfig.DUTPrs(1).RangeMain.GetParent.Extra3  
  
' Read A/F C0 from PMM; make unit conversions, write to calcoef1:
C0raw = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd1), False)
cDebug.LogStatus "the raw C0 string is: " & C0raw
C0 = qextract(C0raw,1,2,",")
    If ModelTypeID = 5 Then 'this is a PM500, therefore the C0 is in units of kPa
        cdebug.LogStatus "the A/F C0 value in kPa is: " & C0
        C0cnvt = cCOMPASS.UnitConversion(CDbl(C0),CInt(DUTUnit),3,0)
        cdebug.LogStatus "the A/F C0 value converted to the DUT unit of measure is: " &
C0cnvt & DUTunitText
    Else
        'PM200 and PM600 always read the C0 in units of psi
        cDebug.LogStatus "The C0 value before conversion is: " & C0
        C0cnvt = cCOMPASS.UnitConversion(CDbl(C0),CInt(DUTUnit),9,0)
        cdebug.LogStatus "the A/F C0 value converted to the DUT unit of measure is: " &
C0cnvt & DUTunitText
    End If
    cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient1 = C0cnvt  
  
' Read A/F C1 from DUT; write to calcoef2
C1raw = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd2), False)
C1 = qextract(C1raw,1,2,",")
cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient2 = C1
cDebug.LogStatus "As Found C1 coefficient: " & C1  
  
' Read A/F zOffset from PMM, write to calcoef3:
zOffset = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd3), False)
cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient3 = zOffset
cDebug.LogStatus "As Found zOffset is: " & zOffset
cDebug.LogStatus "the zOffset value is in units of: " &
cConfig.DUTPrs(1).RangeMain.UnitFinalText  
  
' Read the last saved cal edit date:
calDate = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd4), False)
cCOMPASS.DataCollection(1).DUT.LastCalDate = calDate
cDebug.LogStatus "As Found last cal date: " & calDate  
  
' Place the 6270A into MEASURE mode to prepare for applying the pressure:
cConfig.DUTPrs(1).IoSendCommand "OUTP:MODE MEASURE",True,0
TimeDelay 2
mode = cConfig.DUTPrs(1).IoSendCommand("OUTP:MODE?",False,0)
```

```

        cdebug.LogStatus "Measurement mode: " & mode
        cdebug.LogStatus "Placed the 6270A into MEASURE mode"

        cdebug.LogStatus "<End of CurrentTestStep 1010>"

'<The test definition will now run the first cycle of points>
'<Continue with pressure cycle 1 (as-found data)>

'Vent the 6270A and allow it to run AutoZero...
ElseIf cCOMPASS.CurrentTestStep = 1100 Then
    cDebug.LogStatus "starting step <1100>..."
    DAQ = cCOMPASS.cConfig.SetPrs(1).RangeMain.GetParent.DaqType

    If DAQ <> 0 Then
        'Vent both the Reference and the chassis (DUT). The chassis automatically zeros the
PMM when vented:
        cCOMPASS.cConfig.SetPrs(1).IoSetOutput 0,0,1
        cCOMPASS.cConfig.DUTPrs(1).IoSetOutput 0,0,1
        cDebug.LogStatus "Sent ioSetOutput 0,0,1 to vent the Ref Device."
        cDebug.LogStatus "Sent ioSetOutput 0,0,1 to vent the chassis (DUT)."
        TimeDelay 3
        cCOMPASS.StatusDisplay "Waiting for AutoZero of the PMM..."
        TimeDelay 15
    ElseIf DAQ = 0 Then
        'DaqType mode = 0: manual Set device
        msg = MsgBox("Press the VENT button on the chassis front panel and vent the reference
pressure device. Press OK when fully vented.",0+64+4096,"Manual Vent")
        TimeDelay 3
    End If

    ' Place the chassis into MEASURE mode to prepare for applying the pressure:
    cConfig.DUTPrs(1).IoSendCommand "OUTP:MODE MEASURE",True,0
    TimeDelay 2
    mode = cConfig.DUTPrs(1).IoSendCommand("OUTP:MODE?",False,0)
    cdebug.LogStatus "Placed the chassis into MEASURE mode"

    cDebug.LogStatus "End of step <1100>."

'<currentteststep 1150 - pressure cycle 1 is done, run polyfit>
ElseIf cCOMPASS.CurrentTestStep = 1150 And iC = 1 Then
    cCOMPASS.StatusDisplay "Calculating Adjustment..."
    cDebug.LogStatus "First Cycle is complete. Starting <step 1150> - run PolyFit"
    cDebug.LogStatus "#Data files: " & cCOMPASS.DataCollection.Count

    For i = 1 To cConfig.DUTPrs.Count
        pressurePoints = cCOMPASS.DataCollection(i).NumberOfPressurePoints
        cDebug.LogStatus "DUT: " & i & " - Pressure Points: " & pressurePoints

        Redim Reference_Data(pressurePoints-1)
        Redim DUT_Raw_Data(pressurePoints-1)
        Redim C1_Corrected_Data(pressurePoints-1)

        'populate PMM A/F coefficients from memory:
        C0 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient1
        C1 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient2
        zOffset = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient3
        cDebug.LogStatus "As Found C0 coef: " & C0
        cDebug.LogStatus "As Found C1 coef: " & C1
        cDebug.LogStatus "As Found zOffset: " & zOffset

        ' Iterate through each pressure point
        cDebug.LogStatus "Calculating DUT_Raw(factory)_Data for each test point..."
        For j = 1 To pressurePoints
            ix = CInt(j)-1

            DUTPressure = cCOMPASS.DataCollection(i).DataPointRef(1, 1, 1, CInt(j)).DUTPressure
            Reference_Data(ix) = cCOMPASS.DataCollection(i).DataPointRef(1, 1, 1,
CInt(j)).RefPressure

            'Back out the As Found coefficients from each pressure point to determine the "Factory DUT
Pressure"...
            DUT_Raw_Data(ix)= (DUTPressure-C0-zOffset)/C1
            cDebug.LogStatus "DUT_Raw_Data: " & DUT_Raw_Data(ix)

        Next

        Call Poly_Fit(Reference_Data, DUT_Raw_Data, res, 1)
        slope = res(0,1)
        cDebug.LogStatus "A/L slope: " & slope
        If slope = 0 Then
            new_C1 = 0
        Else

```

```

        new_C1 = 1 / slope
    End If
    cDebug.LogStatus "The new A/L C1 coef " & new_C1

    'Calculate the C1_corrected for each point:
    For j = 0 To pressurePoints - 1
        C1_Corrected_Data(CInt(j)) = DUT_Raw_Data(CInt(j)) * new_C1
        'cDebug.LogStatus "A/L Adjust C1_Corrected_Data(" & j & "): " & C1_Corrected_Data(CInt(j))
    Next

    Call Poly_Fit(Reference_Data, C1_Corrected_Data, res, 1)
    new_C0 = 0 - res(0,0)
    cDebug.LogStatus "The new A/L C0 coef in the DUT unit of measure: " & new_C0

    'Need to convert the new_C0 back to kPa if this is a PM500 calibration. If PM200, 600 then convert
    'the new_C0 back to psi...
    ModelTypeID = cConfig.DUTPrs(1).RangeMain.GetParent.Extra3
    DUTunitText = cConfig.DUTPrs(1).RangeMain.UnitFinalText
    DUTunit = cConfig.DUTPrs(1).RangeMain.UnitFinal
    If ModelTypeID = 5 Then 'this is a PM500, therefore the C0 must go back to units of
    kPa.
        cdebug.LogStatus "the Model Type ID is: " & ModelTypeID
        cdebug.LogStatus "the A/L new_C0 value is: " & new_C0 & DUTunitText
        new_C0bkcnvt =
    cCOMPASS.UnitConversion(CDb1(new_C0),3,CInt(DUTunit),0)
        cdebug.LogStatus "the A/L new_C0 value converted back to kPa is: " &
    new_C0bkcnvt
    Else
        'PM200 and PM600 always read the C0 in units of psi
        new_C0bkcnvt =
    cCOMPASS.UnitConversion(CDb1(new_C0),9,CInt(DUTunit),0)
        cDebug.LogStatus "the A/L new_C0 value converted back to psi is: " &
    new_C0bkcnvt
    End If

    'Writing new data to the PMM. Construct the new write commands, then loop through a Do-Loop
    'to make sure that the new values were written.

    'Repopulate the A/F coefficients:
    C0 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient1
    C1 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient2
    zOffset = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient3
    calDate = cCOMPASS.DataCollection(i).DUT.LastCalDate
        cdebug.LogStatus "the As Found C0 coef is: " & C0
        cdebug.LogStatus "the As Found C1 coef is: " & C1

    'Determine new calibration date:
    new_date = ChassisCal_Date
    cDebug.LogStatus "Calibration & Adjustment date: " & new_date

    'Build new write commands:
    PMMSlotAddress = cConfig.DUTPrs(i).RangeMain.GetParent.Extra2
    cmd1 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1?"
    cmd2 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2?"
    cmd3 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
    cmd4 = "CAL:PRES" & PMMSlotAddress & ":DATE?"

    cmd5 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1 " & new_C0bkcnvt
    cmd6 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2 " & new_C1
    cmd7 = "CAL:PRES" & PMMSlotAddress & ":DATE " & new_date
    cmd8 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA 0"
        cdebug.LogStatus "cmd5 is: " & cmd5
        cdebug.LogStatus "cmd6 is: " & cmd6
        cdebug.LogStatus "cmd7 is: " & cmd7
        cdebug.LogStatus "cmd8 is: " & cmd8

    'Write new coefficients to data file:
    cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient4 = new_C0
    'Note: CalCoef4 is written in units of the DUT Test Unit (typically psi)
    cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient5 = new_C1
    cCOMPASS.DataCollection(i).DUT.CalDueDate = CDate(new_date)

    Do
        'Enable the CAL Mode:
        cConfig.DUTPrs(i).IoSendCommand "CAL:MODE 1",True
        TimeDelay 1
        mode = cConfig.DUTPrs(i).IoSendCommand ("CAL:MODE?",False)
        cdebug.LogStatus "Cal Mode is: " & mode

        'Write the new coeffs:
        Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd5),True)
        TimeDelay 0.5
        Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd6),True)
        TimeDelay 0.5

```

```

    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd7),True)
    TimeDelay 0.5
    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd8),True)
    TimeDelay 0.5

    'Read back the written coeffs:
    WtC0raw = cConfig.DUTPrs(i).IoSendCommand (CStr(cmd1),False)
    TimeDelay 0.5
    WtClraw = cConfig.DUTPrs(i).IoSendCommand (CStr(cmd2),False)
    TimeDelay 0.5
    Wtz = cConfig.DUTPrs(1).IoSendCommand (CStr(cmd3),False)
    timedelay 0.5
        cdebug.LogStatus "WtC0raw value is: " & WtC0raw
        cdebug.LogStatus "WtClraw value is: " & WtClraw
        cdebug.LogStatus "Wtz value is: " & Wtz
    'Remove the leading coef identifier to attain only the value...
    WtC0 = qextract(WtC0raw,1,2,"")
    WtCl = qextract(WtClraw,1,2,"")
        cDebug.LogStatus "WtC0 value is: " & WtC0
        cDebug.LogStatus "WtCl value is: " & WtCl
        cDebug.LogStatus "Wtz value is: " & Wtz

    'Check if the A/L is different from the A/F (it should be):
    C0Diff = C0 - WtC0
    C1Diff = C1 - WtCl
    ZoffDiff = zOffset - Wtz
        cDebug.LogStatus "C0 Difference: " & C0Diff
        cDebug.LogStatus "C1 Difference: " & C1Diff
        cDebug.LogStatus "Zoff Difference: " & ZoffDiff

    If cCOMPASS.SystemAbort Then Exit Function
    If C0Diff <> 0 And C1Diff <> 0 Then Exit Do
    cdebug.LogStatus "End of Loop."
    cCOMPASS.StatusDisplay "Writing new coefficients..."
Loop
Next

Set obj = Nothing
cdebug.LogStatus "Finished with CurrentTestStep 1150"

'<currentteststep 2000 - test is complete, A/L data has been collected. Save the coeffs>

ElseIf cCOMPASS.CurrentTestStep = 2000 Then
cDebug.LogStatus "Test complete. CurrentTestStep 2000. Get final zOffset"

'Confirm the final zOffset is zero:
PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2
cmd9 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
zOffset = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd9), False)
cDebug.LogStatus "Verification of A/L zOffset: " & zOffset

' Write the final zOffset value to the data file:
cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient6 = zOffset

*****
'This block of code is used to check the fit of the As Left corrections. The key variable is named "TOLFACT".
'If TOLFACT = 1 then the As Left data must fit inside the 1 year tolerance specification. If set to <1 then
'this acts as a form of guardbanding.

TOLFACT = .70 '70% gaurdband. User definable.

Dim MessageBox1
Dim MessageBox2
For i = 1 To cCOMPASS.DataCollection(1).NumberofPressurePoints
    For j = 1 To cConfig.DUTPrs.Count
        Ref = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).RefPressure
        DUT = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).DUTPressure
        Tol = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).Tolerance
        ST = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).Status

        If abs(DUT-Ref) > Tol * TOLFACT Then
            FAIL = True

        'Give user the option of keeping or restoring the A/F calibration:
        cCOMPASS.StatusDisplay "Waiting for the user make a decision on the As Left
results..."
```

MessageBox1 = msgbox("One or more of the As Left data points fall outside of the allowed tolerance. Click 'Yes' to accept the results, or click 'No' to cancel the adjustment and restore the As Found coefficients.",4+48+4096,"Accept Results?")

If MessageBox1 = 7 Then 'Note: 7 = yes, 6 = no

```

        Messagebox2 = msgbox("This will delete the newly determined As Found
coefficients and return the PMM to how it was prior to the calibration. Click 'Yes' to restore.",4+32+4096,
"Are you sure?")
        cdebug.LogStatus "The user chose to revert back to the A/F coefficients."
        End If
        If Messagebox2 = 6 Then 'Note: 6 = yes, I'm sure
        RestoreCal

        msgbox "The As Found coefficients have been restored."
        End If

        Exit For
        End If
    Next
    If fail = True Then Exit For

'END OF GUARD BANDING CODE.
*****
```

Next

'Send a SAVE command for chassis with old firmware - pre-v1.10 (6270A), v1.02 (8270A)...

ID = cConfig.DUTPrs(1).IoSendCommand("*IDN?", False)

cdebug.LogStatus "Response from ID command is: " & ID

ModelRaw = qextract(ID,1,2,";")

cdebug.LogStatus "The chassis model is: "& ModelRaw

Model = Left(ModelRaw,4)

cdebug.LogStatus "the simplified model nomenclature is: " & Model

Firmware = qextract(ID,3,4,";")

FirmwareCdbl = Cdbl(Firmware)

cdebug.LogStatus "The chassis firmware rev is: " & FirmwareCdbl

PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2

cmd10 = "CAL:PRES" & PMMSlotAddress & ":SAVE"

If Model = "6270" And FirmwareCdbl < 1.10 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 6270A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

ElseIf Model = "8270" And FirmwareCdbl < 1.02 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 8270A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

ElseIf Model = "8370" And FirmwareCdbl < 1.02 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 8370A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

End If

End If

End Function

ChassisCal_Absolute:

```

'*****
'Used with ChassisCal_Master to calibrate and adjust Absolute mode PMMs. Applicable for 6270A / 8270A / 8370A.
'Called by the master test macro "ChassisCal_Master". The Absolute mode calibrations do not zero the PMM prior
'To taking A/F Data.

'*****

Function ChassisCal_Absolute(iT, iL, iC, iP, cTest, cConfig)
Dim DUT_Raw_Data()
Dim Reference_Data()
Dim C1_Corrected_Data()
Dim res(4,7)

C0 = 0
C1 = 0
zOffset = 0

'<currentteststep 1010 - record AF coeffs.>
If cCOMPASS.CurrentTestStep = 1010 Then ' Test Definition data files were created for each DUT.
    cDebug.LogStatus "Starting <step 1010>"

    'Ensure the Chassis is in Absolute mode:
    cConfig.DUTPrs(1).IoSendCommand "SENS:MODE ABS",True

    DAQ = cCOMPASS.cConfig.SetPrs(1).RangeMain.GetParent.DaqType
    cDebug.LogStatus "Data Acquisition Type for SetPrs is: " & DAQ

    If DAQ <> 0 Then
        'Vent both the Reference and 6270A.
        cCOMPASS.cConfig.SetPrs(1).IoSetOutput 0,0,1
        cCOMPASS.cConfig.DUTPrs(1).IoSetOutput 0,0,1
            cDebug.LogStatus "Sent ioSetOutput 0,0,1 to vent the Ref Device."
            cDebug.LogStatus "Sent ioSetOutput 0,0,1 to vent the 6270A."
        cCOMPASS.StatusDisplay "Venting the pressure..."
        TimeDelay 3
    ElseIf DAQ = 0 Then
        'DagType mode = 0: manual Set device
        msg = MsgBox("Press the VENT button on the 6270A front panel and vent the reference pressure
device. Press OK when fully vented.",0+64+4096,"Manual Vent")
        TimeDelay 5
    End If

    cDebug.LogStatus "Capture As Found Coeffs..."
    PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2
    cDebug.LogStatus "The PMM slot address is: " & PMMSlotAddress

    cmd1 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1?"
    cmd2 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2?"
    cmd3 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
    cmd4 = "CAL:PRES" & PMMSlotAddress & ":DATE?"

    'Identify the unit of measure being used for the test and the PMM Type:
    DUTunit = cConfig.DUTPrs(1).RangeMain.UnitFinal
    DUTunitText = cConfig.DUTPrs(1).RangeMain.UnitFinalText
    ModelTypeID = cConfig.DUTPrs(1).RangeMain.GetParent.Extra3

    ' Read A/F C0 from PMM; make unit conversions, write to calcoef1:
    C0raw = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd1), False)
    cDebug.LogStatus "the raw C0 string is: " & C0raw
    C0 = qextract(C0raw,1,2,",")
        If ModelTypeID = 5 Then 'this is a PM500, therefore the C0 is read in units of kPa
            cDebug.LogStatus "the A/F C0 value in kPa is: " & C0
            C0cnvt = cCOMPASS.UnitConversion(CDbl(C0),CInt(DUTunit),3,0)
                cDebug.LogStatus "the A/F C0 value converted to the DUT unit of measure is: " &
C0cnvt & DUTunitText
        Else
            'PM200 and PM600 always read the C0 in units of psi
            cDebug.LogStatus "The C0 value before conversion is: " & C0
            C0cnvt = cCOMPASS.UnitConversion(CDbl(C0),CInt(DUTunit),9,0)
                cDebug.LogStatus "the A/F C0 value converted to the DUT unit of measure is: " &
C0cnvt & DUTunitText
        End If
    cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient1 = C0cnvt

    ' Read A/F C1 from DUT; write to calcoef2
    C1raw = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd2), False)
    C1 = qextract(C1raw,1,2,",")
    cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient2 = C1
        cDebug.LogStatus "As Found C1 coefficient: " & C1

```

```

' Read A/F zOffset from PMM, write to calcoef3:
' Note: zOffset is displayed on front panel and queried in the same unit of measure as the chassis.
zOffset = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd3), False)
cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient3 = zOffset
    cDebug.LogStatus "As Found zoffset is: " & zOffset
    cDebug.LogStatus "the zOffset value is in units of: " &
cConfig.DUTPrs(1).RangeMain.UnitFinalText

' Read the last saved cal edit date:
calDate = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd4), False)
cCOMPASS.DataCollection(1).DUT.LastCalDate = calDate
    cDebug.LogStatus "As Found last cal date: " & calDate

' Place the 6270A into MEASURE mode to prepare for applying the pressure:
cConfig.DUTPrs(1).IoSendCommand "OUTP:MODE MEASURE",True,0
TimeDelay 2
mode = cConfig.DUTPrs(1).IoSendCommand("OUTP:MODE?",False,0)
    cdebug.LogStatus "Measurement mode: " & mode
    cdebug.LogStatus "Placed the 6270A into MEASURE mode"

cdebug.LogStatus "<End of CurrentTestStep 1010>"

'<The test definition will now run the first cycle of points>
'<Continue with pressure cycle 1 (as-found data)>

'<currentteststep 1150 - pressure cycle 1 is done, run polyfit>
ElseIf cCOMPASS.CurrentTestStep = 1150 And iC = 1 Then
cCOMPASS.StatusDisplay "Calculating Adjustment..."
    cDebug.LogStatus "First Cycle is complete. Starting <step 1150> - run PolyFit"
    cDebug.LogStatus "#Data files: " & cCOMPASS.DataCollection.Count

For i = 1 To cConfig.DUTPrs.Count
    pressurePoints = cCOMPASS.DataCollection(i).NumberofPressurePoints
    cDebug.LogStatus "DUT: " & i & " - Pressure Points: " & pressurePoints

    Redim Reference_Data(pressurePoints-1)
    Redim DUT_Raw_Data(pressurePoints-1)
    Redim C1_Corrected_Data(pressurePoints-1)

    'populate PMM A/F coefficients from memory:
    C0 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient1
    C1 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient2
    zOffset = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient3
        cDebug.LogStatus "As Found C0 coef: " & C0
        cDebug.LogStatus "As Found C1 coef: " & C1
        cDebug.LogStatus "As Found zOffset: " & zOffset

    'Iterate through each pressure point:
        cDebug.LogStatus "Calculating DUT_Raw(factory)_Data for each test point..."
    For j = 1 To pressurePoints
        ix = CInt(j)-1

        DUTPressure = cCOMPASS.DataCollection(i).DataPointRef(1, 1, 1, CInt(j)).DUTPressure
            Reference_Data(ix) = cCOMPASS.DataCollection(i).DataPointRef(1, 1, 1,
CInt(j)).RefPressure

        'Back out the As Found coefficients from each pressure point to determine the "Factory DUT
        Pressure"...
            DUT_Raw_Data(ix)= (DUTPressure-C0-zOffset)/C1
            cDebug.LogStatus "DUT_Raw_Data: " & DUT_Raw_Data(ix)

        Next

        Call Poly_Fit(Reference_Data, DUT_Raw_Data, res, 1)
        slope= res(0,1)
            cDebug.LogStatus "A/L slope: " & slope
        If slope = 0 Then
            new_C1 = 0
        Else
            new_C1 = 1 / slope
        End If
            cDebug.LogStatus "The new A/L C1 coef " & new_C1

        'Calculate the C1_corrected for each point:
        For j = 0 To pressurePoints - 1
            C1_Corrected_Data(CInt(j)) = DUT_Raw_Data(CInt(j)) * new_C1
            'cDebug.LogStatus "A/L Adjust C1_Corrected_Data(" & j & "): " & C1_Corrected_Data(CInt(j))
        Next

        Call Poly_Fit(Reference_Data, C1_Corrected_Data, res, 1)
        new_C0 = 0 - res(0,0)

```

```

cDebug.LogStatus "The new A/L C0 coef in the DUT unit of measure: " & new_C0

'Need to convert the new_C0 back to kPa if this is a PM500 calibration. If PM200, 600 then convert
the new_C0 back to psi...
DUTunit = cConfig.DUTPrs(1).RangeMain.UnitFinal
DUTunitText = cConfig.DUTPrs(1).RangeMain.UnitFinalText
ModelTypeID = cConfig.DUTPrs(1).RangeMain.GetParent.Extra3

If ModelTypeID = 5 Then 'this is a PM500, therefore the new_C0 must go back to units
of kPa to be written to chassis.
    cdebug.LogStatus "the Model Type ID is: " & ModelTypeID
    cdebug.LogStatus "the A/L new_C0 value is: " & new_C0 & DUTunitText
    new_C0bkcnvt =
cCOMPASS.UnitConversion(Cdbl(new_C0),3,CInt(DUTunit),0)
    cdebug.LogStatus "the A/L new_C0 value converted back to kPa is: " &
new_C0bkcnvt
Else
    'PM200 and PM600 always read the C0 in units of psi
    new_C0bkcnvt =
cCOMPASS.UnitConversion(Cdbl(new_C0),9,CInt(DUTunit),0)
    cDebug.LogStatus "the A/L new_C0 value converted to psi is: " &
new_C0bkcnvt
End If

'Writing new data to the PMM. Construct the new write commands, then loop through a Do-Loop
'to make sure that the new values were written...

'Repopulate the A/F coefficients:
C0 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient1
C1 = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient2
zOffset = cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient3
calDate = cCOMPASS.DataCollection(i).DUT.LastCalDate
    cdebug.LogStatus "the As Found C0 coef is: " & C0
    cdebug.LogStatus "the As Found C1 coef is: " & C1

'Determine new calibration date:
new_date = ChassisCal_Date
cDebug.LogStatus "Calibration & Adjustment date: " & new_date

'Build new write commands:
PMMSlotAddress = cConfig.DUTPrs(i).RangeMain.GetParent.Extra2
cmd1 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1?"
cmd2 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2?"
cmd3 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
cmd4 = "CAL:PRES" & PMMSlotAddress & ":DATE?"

cmd5 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL1 " & new_C0bkcnvt
cmd6 = "CAL:PRES" & PMMSlotAddress & ":DATA:VAL2 " & new_C1
cmd7 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA 0"
cmd8 = "CAL:PRES" & PMMSlotAddress & ":DATE " & new_date

    cdebug.LogStatus "cmd5 is: " & cmd5
    cdebug.LogStatus "cmd6 is: " & cmd6
    cdebug.LogStatus "cmd7 is: " & cmd7
    cdebug.LogStatus "cmd8 is: " & cmd8

'Write new coefficients to data file:
cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient4 = new_C0
'Note: CalCoef4 is written in units of the DUT Test Unit (typically psi)
cCOMPASS.DataCollection(i).DUT.CalibrationCoefficient5 = new_C1
cCOMPASS.DataCollection(i).DUT.CalDueDate = CDate(new_date)

Do
    'Enable the CAL Mode:
    cConfig.DUTPrs(i).IoSendCommand "CAL:MODE 1",True
    TimeDelay 1
    mode = cConfig.DUTPrs(i).IoSendCommand ("CAL:MODE?",False)
    cdebug.LogStatus "Cal Mode is: " & mode

    'Write the new coefs:
    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd5),True)
    TimeDelay 0.5
    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd6),True)
    TimeDelay 0.5
    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd7),True)
    TimeDelay 0.5
    Call cConfig.DUTPrs(i).IoSendCommand(CStr(cmd8),True)
    TimeDelay 0.5

    'Read back the written coeffs:
    WtC0raw = cConfig.DUTPrs(i).IoSendCommand (CStr(cmd1),False)
    TimeDelay 0.5
    WtC1raw = cConfig.DUTPrs(i).IoSendCommand (CStr(cmd2),False)
    TimeDelay 0.5

```

```

WtZ = cConfig.DUTPrs(1).IoSendCommand (CStr(cmd3),False)
TimeDelay 0.5
    cdebug.LogStatus "WtC0raw value is: " & WtC0raw
    cdebug.LogStatus "WtClraw value is: " & WtClraw
    cdebug.LogStatus "Wtz value is: " & Wtz
'Remove the leading coef identifier to attain only the value...
WtC0 = qextract(WtC0raw,1,2,"")
WtCl = qextract(WtClraw,1,2,"")
    cDebug.LogStatus "WtC0 value is: " & WtC0
    cDebug.LogStatus "WtCl value is: " & WtCl
    cDebug.LogStatus "Wtz value is: " & Wtz

'Check if the A/L is different from the A/F (it should be):
C0Diff = C0 - WtC0
C1Diff = C1 - WtCl
ZoffDiff = zOffset - Wtz
    cDebug.LogStatus "C0 Difference: " & C0Diff
    cDebug.LogStatus "C1 Difference: " & C1Diff
    cDebug.LogStatus "Zoff Difference: " & ZoffDiff

If cCOMPASS.SystemAbort Then Exit Function
If C0Diff <> 0 And C1Diff <> 0 Then Exit Do
cdebug.LogStatus "End of Loop."

cCOMPASS.StatusDisplay "Writing new coefficients..."
```

Loop

Next

```

Set obj = Nothing
cdebug.LogStatus "Finished with CurrentTestStep 1150"
```

'<currentteststep 2000 - test is complete, AL data has been collected. Save the coefs>

```

ElseIf cCOMPASS.CurrentTestStep = 2000 Then
cDebug.LogStatus "Test complete. CurrentTestStep 2000. Get final zOffset"

'Confirm the final zOffset is zero:
PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2
cmd9 = "CAL:PRES" & PMMSlotAddress & ":ZERO:DATA?"
zOffset = cConfig.DUTPrs(1).IoSendCommand(CStr(cmd9), False)
cDebug.LogStatus "Verification of A/L zOffset: " & zOffset

' Write the final zOffset value to the data file:
cCOMPASS.DataCollection(1).DUT.CalibrationCoefficient6 = zOffset
```

```

*****  

'This block of code is used to check the fit of the As Left corrections. The key variable is named "TOLFACT".  

'If TOLFACT = 1 then the As Left data must fit inside the 1 year tolerance specification. If set to <1 then  

'this acts as a form of guardbanding.
```

```

TOLFACT = .70      '70% gaurdband. User definable.
```

```

Dim MessageBox1
Dim MessageBox2

For i = 1 To cCOMPASS.DataCollection(1).NumberofPressurePoints
    For j = 1 To cConfig.DUTPrs.Count
        Ref = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).RefPressure
        DUT = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).DUTPressure
        Tol = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).Tolerance
        ST = cCOMPASS.DataCollection(j).DataPointRef(1,1,2,CInt(i)).Status

        If abs(DUT-Ref) > Tol * TOLFACT Then
            FAIL = True

            'Give user the option of keeping or restoring the A/F calibration:
            cCOMPASS.StatusDisplay "Waiting for the user make a decision on the As Left
results..."
```

```

            MessageBox1 = msgbox("One or more of the As Left data points fall outside of
the allowed tolerance. Click 'Yes' to accept the results, or click 'No' to cancel the adjustment and restore
the As Found coefficients.",4+48+4096,"Accept Results?")
```

```

            If MessageBox1 = 7 Then 'Note: 7 = yes, 6 = no
                MessageBox2 = msgbox("This will delete the newly determined As Found
coefficients and return the PMM to how it was prior to the calibration. Click 'Yes' to restore.",4+32+4096,
"Are you sure?")
                cdebug.LogStatus "The user chose to revert back to the A/F coefficients."
            End If
            If MessageBox2 = 6 Then 'Note: 6 = yes, I'm sure
                RestoreCal
            End If
```

```

            msgbox "The As Found coefficients have been restored."
        End If
```

```

        Exit For
    End If
    Next
    If fail = True Then Exit For

'END OF GUARD BANDING CODE.
*****
```

Next

'Send a *SAVE* command for chassis with old firmware - pre-v1.10 (6270A), v1.02 (8270A)...

ID = cConfig.DUTPrs(1).IoSendCommand("*IDN?", False)

cdebug.LogStatus "Response from ID command is: " & ID

ModelRaw = **gextract**(ID,1,2,"")

cdebug.LogStatus "The chassis model is: " & ModelRaw

Model = Left(ModelRaw,4)

cdebug.LogStatus "the simplified model nomenclature is: " & Model

Firmware = **gextract**(ID,3,4,"")

FirmwareCdbl = CDbl(Firmware)

cdebug.LogStatus "The chassis firmware rev is: " & FirmwareCdbl

PMMSlotAddress = cConfig.DUTPrs(1).RangeMain.GetParent.Extra2

cmd10 = "CAL:PRES" & PMMSlotAddress & ":SAVE"

If Model = "6270" And FirmwareCdbl < 1.10 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 6270A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

ElseIf Model = "8270" And FirmwareCdbl < 1.02 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 8270A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

ElseIf Model = "8370" And FirmwareCdbl < 1.02 Then

Call cConfig.DUTPrs(1).IoSendCommand(CStr(cmd10),True)

TimeDelay 0.5

cdebug.LogStatus "the detected model was 8370A with an older ver of firmware."

cdebug.LogStatus "Data saved to PMM."

End If

End If

End Function