



700HPPM Manifold – TEST port seal (not included in 700HPP seal kits)



Original o-ring is urethane size 108 with these dimensions:
OD = 11.3 mm (0.443 in.)
ID = 6.0 mm (0.237 in.)
Thickness = 2.6 mm (0.103 in.)

Notes: You can likely find this at a local o-ring supplier. They know what o-ring size 108 is. Alternately the o-ring or a kit of o-rings might be available through the supplier, FasTest. Their part number is for a kit of 5 is S1108-U70.

We do not offer this o-ring in urethane material but do in Viton.

Fluke Part Number: 3134871

Description: O-RING,FKM,VITON,.103IN W,.237IN ID,75 DURO,BROWN,SIZE 2-108